

---

# **phototdt Documentation**

***Release 0.0.1***

**Matias Andina**

**Apr 18, 2023**



**CONTENTS:**

<b>1</b>	<b>phototdt</b>	<b>1</b>
1.1	Features . . . . .	1
1.2	Credits . . . . .	1
<b>2</b>	<b>Installation</b>	<b>3</b>
2.1	Stable release . . . . .	3
2.2	From sources . . . . .	3
<b>3</b>	<b>Usage</b>	<b>5</b>
<b>4</b>	<b>Contributing</b>	<b>7</b>
4.1	Types of Contributions . . . . .	7
4.2	Get Started! . . . . .	8
4.3	Pull Request Guidelines . . . . .	9
4.4	Tips . . . . .	9
4.5	Deploying . . . . .	9
<b>5</b>	<b>Credits</b>	<b>11</b>
5.1	Development Lead . . . . .	11
5.2	Contributors . . . . .	11
<b>6</b>	<b>History</b>	<b>13</b>
6.1	0.0.1 (2022-10-05) . . . . .	13
<b>7</b>	<b>Indices and tables</b>	<b>15</b>



## **PHOTOTDT**

This Python package contains functions to get photometry data from a Tucker-Davis Technology (TDT) photometry system and calculate dFF using methods developed by Martianova and colleagues. For more information on the analysis method, you can visit [Martianova, E., Aronson, S., Proulx, C.D. Multi-Fiber Photometry to Record Neural Activity in Freely Moving Animal. J. Vis. Exp. \(152\), e60278, doi: 10.3791/60278 \(2019\).](#) Implementation details and other language implementations (R, Matlab) are archived in the [publication's repository](#)

- Free software: BSD license
- Documentation: <https://phototdt.readthedocs.io>.

### **1.1 Features**

This package reads TDT data from the directory of the block (e.g., `photometry_dir`)

- Use `photo_data = phototdt.get_tdt_data(photometry_dir)` to read and obtain a DataFrame photometry data.
- Use `phototdt.tdt_to_csv.tdt_to_csv(photometry_dir)` to convert block to a csv file and calculate zdFF on the 465 channel.
- Use `phototdt.get_cam_timestamps(photometry_dir)` to read camera timestamps from block.

### **1.2 Credits**

This package was created with [Cookiecutter](#) and the [audreyr/cookiecutter-pypackage](#) project template.



## INSTALLATION

### 2.1 Stable release

To install phototdt, run this command in your terminal:

```
$ pip install phototdt
```

This is the preferred method to install phototdt, as it will always install the most recent stable release.

If you don't have [pip](#) installed, this [Python installation guide](#) can guide you through the process.

### 2.2 From sources

The sources for phototdt can be downloaded from the [Github repo](#).

You can either clone the public repository:

```
$ git clone git://github.com/matiasandina/phototdt
```

Or download the [tarball](#):

```
$ curl -OJL https://github.com/matiasandina/phototdt/tarball/master
```

Once you have a copy of the source, you can install it with:

```
$ python setup.py install
```





## USAGE

To use phototdt in a project:

```
import phototdt
session_folder = "path/to/block/folder"
# Get tdt data
photo_data = phototdt.get_tdt_data(session_folder)
# Get camera timestamps
from phototdt.phototdt import get_cam_timestamps
cam_timestamps = get_cam_timestamps(folder=session_folder)
# Convert photometry data from block to csv (interactive if folder is None)
from phototdt.tdt_to_csv import tdt_to_csv
tdt_to_csv(session_folder)
```

::

**You can rename the block contents into BIDS format::**

```
from phototdt.rename_block import rename_block
session_folder = "path/to/block/folder"
rename_block(session_folder)
```

::



## CONTRIBUTING

Contributions are welcome, and they are greatly appreciated! Every little bit helps, and credit will always be given. You can contribute in many ways:

### 4.1 Types of Contributions

#### 4.1.1 Report Bugs

Report bugs at <https://github.com/matiasandina/phototdt/issues>.

If you are reporting a bug, please include:

- Your operating system name and version.
- Any details about your local setup that might be helpful in troubleshooting.
- Detailed steps to reproduce the bug.

#### 4.1.2 Fix Bugs

Look through the GitHub issues for bugs. Anything tagged with “bug” and “help wanted” is open to whoever wants to implement it.

#### 4.1.3 Implement Features

Look through the GitHub issues for features. Anything tagged with “enhancement” and “help wanted” is open to whoever wants to implement it.

#### 4.1.4 Write Documentation

phototdt could always use more documentation, whether as part of the official phototdt docs, in docstrings, or even on the web in blog posts, articles, and such.

### 4.1.5 Submit Feedback

The best way to send feedback is to file an issue at <https://github.com/matiasandina/phototdt/issues>.

If you are proposing a feature:

- Explain in detail how it would work.
- Keep the scope as narrow as possible, to make it easier to implement.
- Remember that this is a volunteer-driven project, and that contributions are welcome :)

## 4.2 Get Started!

Ready to contribute? Here's how to set up *phototdt* for local development.

1. Fork the *phototdt* repo on GitHub.
2. Clone your fork locally:

```
$ git clone git@github.com:your_name_here/phototdt.git
```

3. Install your local copy into a virtualenv. Assuming you have virtualenvwrapper installed, this is how you set up your fork for local development:

```
$ mkvirtualenv phototdt
$ cd phototdt/
$ python setup.py develop
```

4. Create a branch for local development:

```
$ git checkout -b name-of-your-bugfix-or-feature
```

Now you can make your changes locally.

5. When you're done making changes, check that your changes pass flake8 and the tests, including testing other Python versions with tox:

```
$ flake8 phototdt tests
$ python setup.py test or pytest
$ tox
```

To get flake8 and tox, just pip install them into your virtualenv.

6. Commit your changes and push your branch to GitHub:

```
$ git add .
$ git commit -m "Your detailed description of your changes."
$ git push origin name-of-your-bugfix-or-feature
```

7. Submit a pull request through the GitHub website.

## 4.3 Pull Request Guidelines

Before you submit a pull request, check that it meets these guidelines:

1. The pull request should include tests.
2. If the pull request adds functionality, the docs should be updated. Put your new functionality into a function with a docstring, and add the feature to the list in README.rst.
3. The pull request should work for Python 3.5, 3.6, 3.7 and 3.8, and for PyPy. Check [https://travis-ci.com/matiasandina/phototdt/pull\\_requests](https://travis-ci.com/matiasandina/phototdt/pull_requests) and make sure that the tests pass for all supported Python versions.

## 4.4 Tips

To run a subset of tests:

```
$ pytest tests.test_phototdt
```

## 4.5 Deploying

A reminder for the maintainers on how to deploy. Make sure all your changes are committed (including an entry in HISTORY.rst). Then run:

```
$ bump2version patch # possible: major / minor / patch
$ git push
$ git push --tags
```

Travis will then deploy to PyPI if tests pass.



## CREDITS

### 5.1 Development Lead

- Matias Andina <[matiasandina@gmail.com](mailto:matiasandina@gmail.com)>

### 5.2 Contributors

This package contains substantial amount of work from the following contributors:

- Ekaterina Martianova <[ekaterina.martianova.1@ulaval.ca](mailto:ekaterina.martianova.1@ulaval.ca)>
- Renato Lombardo <[renato.lombardo@unipa.it](mailto:renato.lombardo@unipa.it)>





## HISTORY

### 6.1 0.0.1 (2022-10-05)

- First github release.



## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`